

DWG. No.

SHEET 1 OF 1

PH.(989) 631-8324 FAX 631-8362 E-MAIL SALES@NEWTECHTM.COM

## GENERAL NO

_ CHECK WHICH APPLY
INTERIOR  ☐ COLD ROLLED STEEL WITH POWDER COATED FINISH ☐ TYPE #304 OR #316 STAINLESS STEEL ☐ CHEMICAL RESISTANT FUME HOOD LINER (CR-900 WHITE)
EXTERIOR  COLD ROLLED STEEL WITH EPOXY POWDER COAT FINISH (COLOR:)  TYPE #304 OR #316 STAINLESS STEEL CEILING ENCLOSURE REQUIRED (TO 1" ABOVE CEILING) SPECIFY FINISHED CEILING HEIGHT
SASH  STANDARD VERTICAL RISING 1/4" SAFETY GLASS  TYPE "A" (HORIZONTAL/VERTICAL) SLIDING  SLIDING HORIZONTAL  RECESSED DOOR BOTTOM TRACK  RAMP STYLE BOTTOM TRACK
ELECTRIC SERVICES INSTALLED  LIGHT SWITCH  FAN SWITCH (SIGNAL LEG ONLY)  110 VOLT GFCI OUTLETS (QUANTITY)  SPECIALTY OUTLETS (QUANTITY VOLTS)  CORD PASS—THRUS (FRONT TO INSIDE) (QUANTITY)  CORD PASS—THRUS (OUTSIDE TO INSIDE) (QUANTITY)  EXPLOSION PROOF FIXTURES
PLUMBING SERVICES INSTALLED  CUPSINK ( SIDE)  AIR (QUANTITY)  NITROGEN (QUANTITY)  VACUUM (QUANTITY)  COLD WATER (QUANTITY)  DI WATER (QUANTITY)  STEAM (QUANTITY)  EXTERNALLY MOUNTED VACUUM BREAKER (QUANTITY)
LIGHTING (FOR THE INTERIOR OF THE HOOD)  STANDARD TWO BUB FLUORESCENT  OTHER
FUME HOOD BLOWER  NOT REQUIRED  INCLUDED — SEE FAN SPEC SHEET